



Quality Medical Products

GRIP-LOK™ H Hydrocolloid Universal Catheter and Tubing Securement



Designed specially for use with pediatric & neonatal patients with fragile skin, the Grip-Lok H™ hydrocolloid securement device works universally with a large variety of lines and tubes. It is designed to hold patient lines in place, prevent unnecessary pulls and improve patient care – while reducing hospital costs.

Grip-Lok H™ is applied to the patients body wherever a tube or line is to be secured. It is a Gentle, Safe, Simple, and Secure method of keeping valuable lines, tubes, and catheters locked in place.

Gentle: Grip-Lok H™ utilizes a hydrocolloid adhesive that is designed to be gentle on the skin of pediatric and neonatal patients with sensitive skin.

This is extremely important with newborns. Because of the immaturity of the skin, newborns are at higher risk for injury and epidermal stripping from adhesive removal. Epidermal stripping can lead to infection, which then leads to longer hospital stays.

As per the NANN (National Association of Neonatal Nurses) Skin Care Guidelines, "a pectin-based adhesive such as hydrocolloid is the preferred adhesive for immature skin."

Safe: Grip-Lok H™ is latex free and packaged sterile, allowing you to use Grip-Lok H™ close to the insertion site or right in the sterile field.

Simple: Grip-Lok H™ is applied in the same simple manner as a bandage and can be easily placed and adjusted by gloved hands.

Secure: Grip-Lok H™ employs specialized adhesives and a selfgripping fastener to



Keep lines, tubes and catheters locked securely in place.



Grip-Lok H™ Advantages

- Strong grip holds important lines in place, saving time, money and patient comfort
- Latex free
- Easily applied and adjusted, even with active children
- Sterility allows use close to the wound or in the sterile field
- Hydrocolloid adhesive helps prevent damage to the fragile neonatal skin.*
- Helps prevent dislodgement and/or migration, thus reducing the risk of lost lines, duplication of work, and patient discomfort or distress.

*Evidence-based Clinical Practice Guideline: Neonatal Skin Care. Washington, DC: Association of Women's Health, Obstetric and Neonatal Nurses; 2001:25 (#3).

Grip-Lok™ H Universal Applications

Grip-Lok H™ can be universally used with many different types and brands of lines, tubes and catheters ranging in size from 1.0 mm (3/64") to 4.5mm (3/16") including the following:

- PICC Lines
- Central Lines
- Feeding Tubes
- IV Lines
- Surgical Drains
- Chest Tubes



[Download Grip-Lok™ H Hydrocolloid Securement Device Brochure](#)

[Download GRIP-LOK™ H Application Sheet](#)



Ordering Information

Product #	Tubing Size Range	Description
3100ESH	1.0mm - 3.0mm (3/64" - 1/8")	Extra Small Grip-Lok H Sterile, Individually Packaged, 100/bx
3200SH	1.5mm - 4.5mm (1/16" - 3/16")	Small Grip-Lok H Sterile, Individually Packaged, 100/bx

All Zefon medical products are sold exclusively through specialty distributors.

Product availability may vary by country/region. Please contact your local Zefon distributor for product, pricing and availability information in your area. A list of distributors can be found in the **Where to Buy** section on this website.

Accessories & Related Products



GRIP-LOK™ Universal Catheter and Tubing Securement



Zefon Arterial Anchor Bandage

IMA
Independent Medical Associates

11733 66th St N Suite 113
Largo, FL 334773
888.548.4462
FAX: 888.548.1462
email: info@I-MA.com
website: www.I-MA.com