



Quality Medical Products

GRIP-LOK™ Universal Catheter and Tubing Securement



Manufactured catheter securement and stabilization devices are now listed as "preferred" in the new *Infusion Nursing Standards of Practice*.

Grip-Lok™ is a stabilization and securement device strong enough to lock catheters, tubes and lines securely in place, yet versatile enough for almost any securement situation.

Safe: Grip-Lok™ is latex free and is packaged in individual sterile units, allowing you to use the Grip-Lok™ close to the wound site or right in the sterile field. The medical grade hypoallergenic tape allows the patient's skin to breathe, thus reducing the risk of rashes and other skin irritations.

Simple: Grip-Lok™ is applied in the same simple manner as a bandage and can be easily placed and adjusted by gloved hands. Available in 3 sizes, each is used and applied in the same identical manner, thus encouraging standardization of important dressings.

Secure: Grip-Lok™ employs specialized adhesives and a self-gripping fastener to keep lines, tubes, and catheters locked securely in place. Once secured, Grip-Lok™ prevents dislodgment and/or migration, thus reducing the risk of lost lines, duplication of work, and patient discomfort or distress.






Standardized: Grip-Lok™ conforms to the new Infusion Nursing Standards of Practice and can eliminate the need for sutures and unscheduled catheter restarts, thus directly reducing needlestick injuries.

Grip-Lok™ can be used with lines, tubes and catheters ranging in size from 1.5 mm (1/16") to 13mm (1/2") including the following:

- Peripherally Inserted Central Catheters (PICC)
- Central Venous Catheters
- IV Lines
- Foley Catheters
- Feeding Tubes
- Surgical Drains
- Chest Tubes
- Bedside Emergency Equipment
- Any place that catheters, tubes or line needs to be secured



Dr. Z Says.....

Grip-Lok offers significant advantages over tape. It is designed to hold significantly more securely and reduce complications. Every line pull or IV restart that is avoided results in better patient care and reduced costs. An example of a comparison using Grip-Lok vs. Tape in an IV Start Kit is available **HERE**.



Ordering Information

Product #	Tubing Size Range	Description
3200S	1.5mm - 4.5mm (1/16" - 3/16")	Small Grip-Lok Sterile, Individually Packaged, 100/bx
3300M	3.0mm - 8.0mm (1/8" - 5/16")	Medium Grip-Lok Sterile, Individually Packaged, 100/bx
3400L	6mm - 13mm (1/4" - 1/2")	Large Grip-Lok Sterile, Individually Packaged, 100/bx

All Zefon medical products are sold exclusively through specialty distributors.

Product availability may vary by country/region. Please contact your local Zefon distributor for product, pricing and availability information in your area. A list of distributors can be found in the **Where to Buy** section on this website.


Accessories & Related Products



GRIP-LOK™ H Hydrocolloid Universal Catheter and Tubing Securement



Zefon Arterial Anchor Bandage



GRIP-LOK™ Universal Foley Catheter Securement



GRIP-LOK® WA Universal PICC Catheter Securement



GRIP-LOK™ NGH Hydrocolloid Nasal Gastric Securement