

# Catheterization Tray

Multiple procedure tray specifically designed for neonatal use

High quality components help improve patient outcomes



MADE  
IN THE  
U.S.A.



**Footprint**  
medical inc.

# Catheterization Tray

## Protection you can count on

Footprint Medical catheterization trays feature high quality components to allow clinicians to provide optimal neonatal care.

Let us show you the Footprint Medical difference... Talk to a representative today!

- Neonatal general procedure tray especially for umbilical vessel catheterization
- Compact bi-level tray design for ease of storage and transport
- High quality disposable stainless steel instrumentation
- Skin Prep, Umbilical Tape, Stopcocks and Suture provided in easy access compartments



\*Umbilical catheters available separately

### Safety Products

Safety Scalpel and sharps protection designed to lock onto and secure contaminated needles.



### Medication Labels

Color coded Medication/Line labels and stopcocks.



### Neonatal Micro Forceps

Both half and full-curved options provided for dilating small vessels.



### Maximum Barrier Drape

Exclusive "Hybrid" drape provides visualization of neonate and prevents contamination.



## CATHETERIZATION TRAY

REORDER NUMBER	QUANTITY
CATHTRAY	5/case

- (1) Adson Forcep
- (2) Iris Micro Forcep 2-3/4"
- (1) Iris Forcep 4"
- (1) Needle Holder
- (1) Safety Scalpel, #11 Blade
- (1) Hemostat – Straight
- (2) Hemostat – Curved
- (1) Scissors, Blunt-Sharp
- (2) 3-Way Stopcock

- (1) 1mL Syringe
- (2) 3mL Syringe
- (1) 5mL Syringe
- (1) Tape Measure
- (1) 5" Probe
- (1) 20ga Needle, 1-1/2"
- (1) Filter Needle, 19GA
- (1) Swabstick, Iodophor (3-pk)
- (1) CSR Wrap

- (3) Drape, Poly-lined
- (1) Umbilical Tape, 12"
- (1) Hybrid Drape w/Clear Window
- (4) Gauze, 4x4
- (4) Gauze, 2x2
- (1) Silk Suture, 4-0
- (2) Sharps Containers, Point-Lok
- (1) Medication Labels

Exclusively Distributed by IMA:  
To request samples or to order,  
please call:  
**(727) 548-4462**  
www.i-ma.com

**IMA**  
Independent Medical Associates

**Footprint**  
medical inc.

www.footprintmed.com

12727 Cimarron Path  
San Antonio, TX 78249  
210.226.2600

PB-CATHTRAY, REV1 20-037