



MaxPlus EMT Cooler - with - Maxify™ Technology



ECONOMICAL



EFFICIENT



SUSTAINABLE



USA-MADE

Product Page:



Qualified cooler for emergency medical transport & air evac operations

Holds 2 (500 mL) whole blood units

From pre-hospital to OR and ER transfusions, trauma packs and MTP, the MaxPlus EMT cooler transports blood products safely, stores efficiently and stands up to repeated in-field and in-hospital use. All with no special handling or bench-time conditioning.

What is Maxify™ Technology?

MAXIFY™ proprietary materials technologies are revolutionizing the shipping of temperature-sensitive biologics, enabling our leading-edge MAXQ no-bench-time packaging solutions to provide unmatched shipping and cost efficiencies, thermal performance, and TOTAL VALUE.

- No-Bench-Time preconditioning required
- Universal (all-season) protection
- Optimized for quick deployment
- Qualified in operational ambient environments
- Proven MaxPlus Protection - eliminates loss

Sales@packmaxq.com | 405-466-5629
Order at: www.PackMaxQ.com



Product Page:



Maximize Savings!

MaxPlus coolers combine unprecedented thermal performance and convenience to decrease lost blood product and increase savings.

SKU #	Capacity	Dimension	Temperature	Validation
EMT10V24	(x2) 500mL units	10.875"x10.5"x9.5"	1 - 6°C	21 hrs standby, 3 in extreme transport



System Components:

- Two PCM5 Cassettes
- One PCM0 Brick

Thermal Performance Plot: Extreme Summer Example

