Polyurethane Enteral Feeding Tube

Designed to provide secure enteral only connections

Open End

Design with precision drilled side holes, enhances fluid flow and decreases potential feeding stagnation and mucosal irritation.

Easy To Identify Hub is clearly marked with "ORAL" and French size to readily identify the tube as an enteral feeding catheter.





Polyurethane Enteral Feeding Tube

Protection you can count on

Footprint Medical enteral feeding tube products offer unique features that allow clinicians more convenience with accuracy in placement, tracking and maintaining the catheter.

Let us show you the Footprint Medical difference... Talk to a representative today!

- Unique polyurethane tube with silicone overlayed hub designed to provide secure enteral only connections
- Labeling supports ISMP recommendations that all distal ports and tubing of access lines clearly indicate what is the use of the port/line
- Bold, indelible markings every centimeter
- Color-coded for patient safety
- Proximal hubs made of inert Polypropylene for excellent chemical resistance and mechanical stability
- Polyurethane's tensile strength allows thin wall construction designed for increased flow rates and decreased clogging with proper maintenance
- Latex and DEHP free
- Integrated silicone closure plug for easy cleaning and reliable sealing.
- Radiopaque stripe for confirmation of tip placement by x-ray
- Smooth, rounded tip surface designed to reduce mucosal irritation



40CM ENTERAL FEEDING TUBES

40CIVI ENTERAL	FEEDING TUBES		
P1FTS4.0-EO	Polyurethane Feeding Tubes	4.0 Fr.	10/case
P1FTS5.0-EO	Polyurethane Feeding Tubes	5.0 Fr.	10/case
P1FTS6.5-EO	Polyurethane Feeding Tubes	6.5 Fr.	10/case
P1FTS8.0-EO	Polyurethane Feeding Tubes	8.0 Fr.	10/case
60CM ENTERAL	FEEDING TUBES		
P1FTM4.0-EO	Polyurethane Feeding Tubes	4.0 Fr.	10/case
P1FTM5.0-EO	Polyurethane Feeding Tubes	5.0 Fr.	10/case
P1FTM6.5-EO	Polyurethane Feeding Tubes	6.5 Fr.	10/case
P1FTM8.0-EO	Polyurethane Feeding Tubes	8.0 Fr.	10/case
90CM ENTERAL	FEEDING TUBES		
P1FTL5.0-EO	Polyurethane Feeding Tubes	5.0 Fr.	10/case
P1FTL6.5-EO	Polyurethane Feeding Tubes	6.5 Fr.	10/case
P1FTL8.0-EO	Polyurethane Feeding Tubes	8.0 Fr.	10/case
P1FTL10.0-EO	Polyurethane Feeding Tubes (110cm)	10.0 Fr.	10/case
Exclusively Distributed by IMA:			

Exclusively Distributed by IMA: To request samples or to order, please call: (727) 548-4462 www.i-ma.com 7301 124th Ave. Largo, FL 33773





5823 Sebastian Place San Antonio, TX 78249 210.226.2600

www.footprintmed.com PB-P1FT REV2 ECO # 21-021