

# **Table of Contents**



Positioning Aids
Developmentally appropriate
positioning solutions that help
support caregivers' efforts in caring
for the most fragile infants



Calming and Soothing
Innovative, evidence-based
solutions designed to help
comfort even the smallest babies



11

Therapy Support
Simple products that support a baby's therapy



Jaundice Management
Proven jaundice management solutions
based on our developmentally friendly
approach to treating newborns

15

# Positioning Aids

# Better Bendy Bumper





These flexible positioning aids mimic the boundaries of the womb to promote containment and flexion, providing a reflex stimulus for extremity extension and recoil.



#### Better Bendy Bumper

The Bendy Bumper comes in two sizes and has unique molding capabilities to maintain their form until reshaped by the caregiver.

	Product #	Legacy #	Description	Quantity per case
0	989805604581	1029244	Better Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	6
2	989805604571	1029243	Better Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	6

## **Better Bendy Bumper Cover**

Two types of covers are available: washable cloth covers for multiple-patient use and disposable covers for single-patient use.

Product #	Legacy #	Description	Quantity per case	
989805606761	92065-B2	Reusable cloth covers for Better Bendy Bumper 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	12	
989805606771	92065-B3	Disposable covers for Better Bendy Bumper 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	50	
989805606701	92065-A2	Reusable cloth covers for Better Bendy Bumper 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	12	
989805606711	92065-A3	Disposable covers for Better Bendy Bumper 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	50	

# SnuggleUp



The SnuggleUp is a soft, cozy nest that helps provide proper positioning for preterm and ill babies. Its padded footroll and soft, adjustable straps allow movement while gently maintaining appropriate positioning.

## Reusable SnuggleUp

SnuggleUps are made from a soft, washable cotton/polyester interlock knit that can be washed and reused.

Product #	Legacy #	Description	Quantity per case
989805606311	91015-X	Extra small for infants < 900 g (< 2 lb)	12
989805606281	91015-A	Small for infants 800 to 1350 g (1.8 to 3 lb)	12
989805606291	91015-B	Medium for infants 1300 to 2250 g (2.8 to 5 lb)	12
989805606301	91015-C	Large for infants > 2250 g (> 5 lb)	12

## Disposable SnuggleUp

Disposable SnuggleUps are particularly useful for admissions, procedures, transport, or when laundry facilities are not available.

Product #.	Legacy #	Description	Quantity per case
989805607421	97013-XDIS	Extra small for infants < 900 g (< 2 lb)	12
989805607391	97013-ADIS	Small for infants 800 to 1350 g (1.8 to 3 lb)	12
989805607401	97013-BDIS	Medium for infants 1300 to 2250 g (2.8 to 5 lb)	12
989805607411	97013-CDIS	Large for infants > 2250 g (> 5 lb)	12

# Frederick T. Frog



Our versatile, frog-shaped beanbags were designed by a NICU nurse to help achieve clinicians' goals for developmentally supportive positioning. Place "Freddy" around the head, neck, or hips, or use it to position an extremity and provide tactile support. Available in two versions: washable/reusable, or wipeable/reusable with launderable or single-use covers.

### Frederick T. Frog – Washable

The washable Freddy Frog can be processed through traditional laundry services.

Product #	Legacy #	Description	Quantity per case
989805603341	97040	Frederick T. Frog, washable	6

## Frederick T. Frog – Wipeable

The wipeable Freddy Frog does not need to be machine washed, just simply wipe it clean. Choose from a cloth/reusable or a disposable cover.

	Product #	Legacy #	Description	Quantity per case
2	989805635141	1084677	Frederick T. Frog, wipeable with reusable covers	6
	989805635151	1084678	Reusable cover	50
	989805635131	1084679	Disposable cover	50

# **Prone Plus**



The Prone Plus provides ventral support to infants with a comfortable, hourglass-shaped pad. They support newborns in the proper prone position by using gravity to facilitate natural rounding of the shoulders, hand-to-mouth coordination, and head-and-trunk alignment.

#### Prone Plus

Provides standard and consistent prone support while helping to maintain a secure, balanced and flexed position.

Product #	Legacy #	Description	Quantity per case	
989805604811	1034167	Prone Plus with knit cover extra small $< 800 \text{ g} (< 1.8 \text{ lb})$	6	
989805604821	1034168	Prone Plus with knit cover small 800 to 1250 g (1.8 to 2.8 lb)	6	
989805604831	1034169	Prone Plus with knit cover medium 1250 to 2000 g (2.8 to 4.4 lb)	6	
989805604891	1034986	Prone Plus with knit cover large > 2000 g (> 4.4 lb)	6	

#### Prone Plus Covers

Soft, washable knit covers are available for each size of the Prone Plus.

Product #	Legacy #	Description	Quantity per case
989805604841	1034170	Washable knit cover (extra small)	6
989805604851	1034171	Washable knit cover (small)	6
989805604861	1034172	Washable knit cover (medium)	6
989805604901	1034987	Washable knit cover (large)	6

# Squishon 2 and 3 and Squishon Mattress







#### Squishon 2 and 3

Provides an exceptionally soft, pliable surface to support and cradle a baby's head and body. The rectangular shape allows for head movement while maintaining a supportive surface. Made from a soft material and filled with a non-toxic gel.

### **Squishon Mattress**

A gel-filled mattress provides a soft surface designed to help relieve pressure points, protect fragile skin and encourage proper positioning. Caregivers can reposition the baby and shift pressure points without having to handle the infant.

#### Squishon 2 and 3

Squishon 2 to support the extremely low birth weight infant. Squishon 3 for older preemies and babies with more resilient skin.

Product #	Legacy #	Description	Quantity per case
989803196361	91033-S	Squishon 2 $6 \times 9 \text{ in } (15 \times 22.5 \text{ cm})^1$	12
989805606461	91033-3	Squishon 3 6 x 9 in (15 x 22.5 cm) <sup>1</sup>	12
989805606481	91033-A	Disposable covers, Squishon 2, 3 & Gel-E Donuts $7 \times 9$ in (18 x 23 cm)	50

#### **Squishon Mattress**

Pliable, gel-filled mattress that cushions the infant and absorbs vibration.

	Product #	Legacy #	Description	Quantity per case
3	989805603301	91017	Squishon Mattress 12 x 22 in (30 x 55 cm) <sup>1</sup>	1
	989805606321	91017-A	Disposable covers, Squishon Mattress 12.5 x 28 in (31.8 x 71.1 cm)	25

# Gel-E Donut



An exceptionally soft skin combined with a unique gel formulation makes Gel-E Donuts soft, yet supportive. They are designed to help relieve pressure caused by prolonged immobility and other conditions where frequent handling may be contraindicated. Reengineered in 2015 to change to a non-porous skin.

Gel-E Donuts

Comfortable & pliable to support and cradle the baby's head.

Product #	Legacy#	Description	Quantity per case	
989803196371	92025-XS	Extra small 7 in diameter (18 cm), fill weight: $325  \mathrm{g}$ (11.5 oz), purple $^1$	12	
989803196381	92025-S	Small 7 in diameter (18 cm), fill weight: 400 g (14 oz), green <sup>1</sup>	12	
989803196391	92025-M	Medium 7 in diameter (18 cm), fill weight: 450 g (16 oz), peach <sup>1</sup>	12	
989805606481	91033-A	Disposable covers, Squishon 2, 3, & Gel-E Donuts	50	

<sup>1.</sup> Has an 18-month shelf life from the date of manufacture.

# Calming & Soothing

## **Sweet-Ease**

This 24% sucrose and purified water solution is a proven way to calm and soothe infants in the NICU, PICU, or well-baby nursery. It's available in a 15 mL cup, 2 mL vial, or 1 mL vial.



#### Sweet-Ease 15 mL cup

This 15 mL cup is aseptically packaged with a peel-off lid that is suitable for dipping a pacifier or for administration via a dropper or syringe.





#### Sweet-Ease 1 & 2 mL vials

Vials simplify clinical workflow by allowing the administration of sucrose directly on to the tip of an infant's tongue or buccal surface.

#### Sweet-Ease 15 mL cup

Aseptically packed cups give clinicians flexibility for administration of sucrose solution.

	Product #	Legacy #	Description	Quantity per case
1	989805603401	99044	Sweet-Ease natural (15 mL cup) <sup>1</sup>	200

#### Sweet-Ease 1 & 2 mL vials

Conveniently packaged vials are great for quick delivery of sucrose solution to the infant.

	Product #	Legacy #	Description	Quantity per case
2	989805637451	1093977	Sweet-Ease natural (2 mL vial) <sup>2</sup>	250
2	00850018561235	N/A	Sweet-Ease natural (1 mL vial) <sup>2</sup>	250

<sup>1.</sup> Has a one-year shelf life from the date of manufacture.

<sup>2.</sup> Has a two-year shelf life from the date of manufacture.

# Snoedel




The Snoedel, originating from a Dutch family tradition, offers parents a means to comfort and bond with their child even when the baby is in the hospital. The pure sheep wool incorporated into the high-quality flannel doll absorbs the parent's scent and slowly releases it, supporting the sensitive olfactory system of the infant.

#### Snoedel

The Snoedel is a flannel bonding doll and positioning aid.

Product #	Legacy #	Description	Quantity per case
989805607221	96029-B	Snoedel (12 unisex print)	12

# **Woombie Med Pods**



#### Med Pod

Woombie Med Pods provide the comfort and security of a swaddle while still allowing for clinical access and procedural support from birth to six months.



## Sleepy Pod

Woombie Sleepy Pods meet AAP Safe Sleep Standards and promote flexion and comfort for full term healthy newborns.

#### Woombie Med Pod

Versatile swaddle designed to be compatible with clinical care in the NICU.

Product #.	Legacy #	Description	Quantity per case
1 731-0359-1	N/A	Woombie Med Pod, Micro Preemie, 1.5-3 lbs	10
731-0359-2	N/A	Woombie Med Pod, Preemie, 3-5 lbs	10
731-0359-3	N/A	Woombie Med Pod, Newborn, 5-13 lbs	10
731-0359-4	N/A	Woombie Med Pod, Infant, 14-19 lbs	10

### Woombie Sleepy Pod

The Sleepy Pod is the top choice for newborn nurseries, L&D units, gift shops, welcome baskets, and NICU graduates who need a safe way to swaddle.

Product #.	Legacy #	Description	Quantity per case
2 731-0363	N/A	Woombie Sleepy Pod, Newborn, 5-13lbs	10

# **Woombie Med Pods**



#### Gastro Pod

The Gastro Pod is a swaddle with a midsection that opens to provide easy clinical access for gastrointestinal or umbilical needs.



#### Trach Pod

The Trach Pod features a deep V-neck for respiratory and cardiac care or for easy access to an oncology patient's port-a-cath or broviac.

#### Woombie Gastro Pod

Swaddle a bigger baby with a g-button or colostomy bag that may need a bit of warmth and containment during their hospital stay.

	Product #.	Legacy #	Description	Quantity per case
1	731-0361-1	N/A	Woombie Gastro Pod, Preemie, 3-5lbs	10
	731-0361-2	N/A	Woombie Gastro Pod, Newborn, 5-13lbs	10

### Woombie Trach Pod

The Trach Pod was designed with trach babies in mind, but can be used for post-operative cardiac care or Oncology patients with a broviac.

	Product #.	Legacy #	Description	Quantity per case
2	731-0362-1	N/A	Woombie Trach Pod, Preemie, 3-5lbs	10
	731-0362-2	N/A	Woombie Trach Pod, Newborn, 5-13lbs	10

# **Therapy Support**

# Wee Sight Transilluminator



The Wee Sight transilluminator helps to locate the tiny veins in a neonate's small limbs for improved insertion of an IV.



### Wee Sight Transilluminator

The Wee Sight uses two AA alkaline batteries (included) and outputs a light wave length of approximately 629 nm.

Product #	Legacy #	Description	Quantity per case
989805603761	1017920	Wee Sight Transilluminator	1

# Sali-Wipes



Sali-Wipes sterile saline (isotonic) wipes come in convenient foil packs. Nursing organizations recommend use of sterile isotonic saline for cleansing, as alcohol wipes can be damaging to fragile infant skin and may be absorbed into the bloodstream.

### Sali-Wipes

Sali-Wipes are used to clean skin after disinfection & prior to invasive procedures.

Product #	Legacy #	Description	Quantity per case
989805607531	98098-013	Sali-Wipes <sup>1</sup> 4 x 6 in (10 x 15.2 cm)	576

1. Has an two-year shelf life from the date of manufacture.

# **Temperature Trend Indicators**



Temperature Trend Indicators provide non-invasive, rapid temperature data for consistent monitoring of a patient's temperature.

### **Temperature Trend Indicators**

Temperature Trend Indicators are accurate to  $\pm$  1°F / 0.5°C, and are ideal for use in a number of care settings.

Product #	Legacy #	Description	Quantity per case
989805623001	LC105D	LC Accu-Bar Temperature Trend Indicator and dispenser box (one degree thermometer) 94°F-105°F (34.5°-40.5°C)	100

# **BBG Nasal Aspirator**



The BBG Nasal Aspirator allows efficient, non-invasive suction of the nasal cavity. It is a safe, easy-to-use bulb-style aspirator that features an over-insertion guard. Two flexible sizes enable caregivers to provide safe suction in a variety of care environments.

### **BBG Nasal Aspirator**

The BBG Nasal Aspirator facilitates gentle suction to remove excess mucus or secretions from a baby's nasal passage.

Product #	Legacy #	Description	Quantity per case	
989805606951	93032-02	BBG Nasal Aspirator 3.25 in (8.25 cm)	50	
989805606961	93032-S	BBG Nasal Aspirator 4.5 in (11.43 cm)	50	

# Jaundice Management

# WeeSpecs





### WeeSpecs phototherapy mask

WeeSpecs' two-piece design is made from a soft, hypoallergenic, light-blocking material that can be adjusted for a custom fit. The inside frame of the goggles seals around the skeletal socket, reducing pressure and allowing the eyes to open and close naturally.

## WeeSpecs Phototherapy Mask

The WeeSpecs offer protection for haby's eyes with a two-piece adjustable design

Product #	Legacy #	Description	Quantity per case	
989805607631	99071-01	WeeSpecs preemie – fits head circumference < 28 cm (< 11 in)	50	
989805607641	99071-02	WeeSpecs small – fits head circumference < 28 to 34 cm (< 11 to 13.4 in)	50	
989805607651	99071-03	WeeSpecs medium – fits head circumference > 34 cm (> 13.4 in)	50	

## WeeSpecs Supreme Phototherapy Mask

Wee Specs Supreme feature an easy-to-adjust strap and snug fit.

	Product #	Legacy #	Description	Quantity per case
2	989805604201	1022894	WeeSpecs Supreme preemie – fits head circumference < 26 cm (< 10.2 in)	50
2	989805604161	1022481	WeeSpecs Supreme small – fits head circumference < 26 to 34 cm (10.2 to 13.4 in)	50
2	989805604211	1022896	WeeSpecs Supreme medium – fits head circumference > 34 cm (> 13.4 in)	50





### WeeSpecs Supreme phototherapy mask

WeeSpecs Supreme incorporate an anatomically contoured goggle design and dark fabric that provide superb light-blocking protection without placing pressure on the eyes. WeeSpecs are held in place with a layer of breathable, slip-resistant foam that reduces movement and are designed to fit newborns to approximately three weeks old.

# Woombie Bili Pod



The Woombie Bili Pod is a swaddle used during phototherapy treatment with nearly 95% of light transmission through the mesh fabric.

### Woombie Bili Pod

The Bili Pod meets the needs of babies undergoing phototherapy, withdrawals (NAS), or cooling therapy.

Product #	Legacy#	Description	Quantity per case	
731-0360-1	N/A	Woombie Bili Pod, Micro Preemie, 1.5-3 lbs	10	
731-0360-2	N/A	Woombie Bili Pod, Preemie, 3-5 lbs	10	
731-0360-3	N/A	Woombie Bili Pod, Newborn, 5-13 lbs	10	

### MRI Compatible Woombie Bili Pod

A Bili Pod designed specifically to keep the infant swaddled during a MRI procedure.

Product #	Legacy #	Description	Quantity per case
731-0567-1	N/A	MRI Compatible Woombie Bili Pod, Micro Preemie, 1.5-3 lbs	10
731-0567-2	N/A	MRI Compatible Woombie Bili Pod, Preemie, 3-5 lbs	10
731-0567-3	N/A	MRI Compatible Woombie Bili Pod, Newborn, 5-13 lbs	10



## For more information visit www.int-bio.com or email sales@int-bio.com.

#### Worldwide

International Biomedical 8206 Cross Park Drive Austin, TX 78754 USA

512.873.0033 sales@int-bio.com

## Europe

Glütschbachstrasse 100 CH-3661 Uetendorf, Switzerland

+41.33.345.66.00 welcome@int-bio.ch



International Biomedical 8206 Cross Park Drive Austin, TX 78754 USA

© 2022 International Biomedical. All rights are reserved.

International Biomedical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Trademarks are the property of International Biomedical or their respective owners.

