



**MaxPlus**  
MTP Cooler®  
-with-  
**Maxify™**  
Technology



ECONOMICAL



EFFICIENT



SUSTAINABLE



USA-MADE

Product Page:



All-in-one Massive Transfusion Protocol Cooler.  
Trusted by hundreds of blood banks globally.

Holds RBCs, thawed plasma, and platelets in one easy to pack cooler. Qualified for up to 12 hours.

This cooler transports blood products safely from OR/ER transfusions, trauma packs, STAT blood, and bedside storage. This dynamic cooler stands up to daily in-hospital use with no special handling or bench-time conditioning.

- One cooler replaces up to three coolers for each MTP request.
- Updated platelet pouch design for optimal thermal performance.
- Highly durable, long-lasting and easy to clean/disinfect.
- Retractable handle and wheels.

Sales@packmaxq.com | 405-466-5629  
Order at: [www.PackMaxQ.com](http://www.PackMaxQ.com)



“MaxPlus is so simple and so easy to use. We’ve decreased our blood product waste drastically since implementation”

- Harbour UCLA Medical Center



Product Page:



### Maximize Savings!

MaxPlus coolers combine unprecedented thermal performance and convenience to decrease lost blood product and increase savings.


SKU #	Capacity	Dimension	Weight
MTP18X12	1-6 RBC + 1-6 PL + 1-2 PLT	18.8" x 18.3" x 15"	21 lbs.



### System includes:

- Three BPOP gel packs
- Two S6 gel packs
- One SG22 gel pack
- Rolling carrier for easy transport

### Ambient Hospital Test - RBC & Warm Plasma



### Ambient Hospital Test - Platelet

