


Umbilical Catheters

Protection you can count on

Footprint Medical umbilical products offer unique features that allow clinicians more convenience with accuracy in placement, tracking and maintaining the catheter.

Let us show you the Footprint Medical difference... Talk to a representative today!

- Made of soft, biocompatible medical grade silicone or polyurethane
- Latex Free and DEHP Free
- Bold, indelible markings every centimeter allow for accurate placement and assessment
- Radiopaque stripe for confirmation of tip placement by x-ray
- Proximal hubs made of inert Polypropylene for excellent chemical resistance and mechanical stability
- Smooth, rounded tip surface designed to reduce thrombogenicity



SINGLE LUMEN UMBILICAL CATHETERS

	P1UVC-2530	Polyurethane Umbilical Catheter	2.5 Fr.	10/case
	P1UVC-3540	Polyurethane Umbilical Catheter	3.5 Fr.	10/case
	P1UVC-4340	Polyurethane Umbilical Catheter	4.3 Fr.	10/case
	P1UVC-540	Polyurethane Umbilical Catheter	5.0 Fr.	10/case
	S1UVC-3540	Silicone Umbilical Catheter	3.5 Fr.	10/case
	S1UVC-540	Silicone Umbilical Catheter	5.0 Fr.	10/case
	S1UVC-6540	Silicone Umbilical Catheter	6.5 Fr.	10/case

DUAL LUMEN UMBILICAL CATHETERS

P2UVC3.5	Polyurethane Umbilical Catheter	3.5 Fr.	10/case
P2UVC5.0	Polyurethane Umbilical Catheter	5.0 Fr.	10/case
S2UVC-436	Silicone Umbilical Catheter	4.0 Fr.	10/case
S2UVC-536	Silicone Umbilical Catheter	5.0 Fr.	10/case

Exclusively Distributed By **IMA**: To request samples or to order, please call:

(727) 548-4462 www.i-ma.com 7301 124th Ave. Largo, FL 33773 USA



